

# MS-PA

## LED Digital Magnetic Stirrer

### Features

- Digital speed control within the range of 0-1500rpm
- Maximum stirring quantity of H<sub>2</sub>O at 3L
- LED display shows stirring speed (MS-PA)



| Specifications                               | MS-PA              | MS-PB              |
|----------------------------------------------|--------------------|--------------------|
| Work plate Dimension                         | Ø135mm (5 inch)    | Ø135mm (5 inch)    |
| Work plate material                          | Nylon+GF           | Nylon+GF           |
| Motor type                                   | Brushless DC motor | Brushless DC motor |
| Stirring positions                           | 1                  | 1                  |
| Max. stirring quantity[H <sub>2</sub> O]     | 3L                 | 3L                 |
| Max.magnetic bar[length]                     | 50mm               | 50mm               |
| Speed range                                  | 200-1500rpm        | 0-1500rpm          |
| Speed display                                | LED                | scale              |
| Protection class                             | IP42               | IP42               |
| Motor output Power                           | 3W                 | 3W                 |
| Power                                        | 15W                | 10W                |
| Voltage,Frequency                            | 100-240V, 50/60Hz  | 100-240V, 50/60Hz  |
| Dimension[WxDxH]                             | 150x260x80mm       | 150x260x80mm       |
| Weight                                       | 1.8kg              | 1.8kg              |
| Permissible ambient temperature and humidity | 10-40°C,80%RH      | 5-40°C,80%RH       |



Strong magnetism



Chemical resistance



LED display